

Neutralizing Rinse

Product Description

Neutralizing rinse is formulated to be used as a neutralizer in the first rinse cycle of instrument cleaning. It assists in keeping instruments and equipment free of mineral build up. It also eliminates alkaline and hard water residue which results in less spotting and prolonged instrument life.

Neutralizing rinse is safe on stainless steel instruments and machines. It is NOT for use on lensed instruments, chrome, or silver instruments.

Applications

- Stainless steel surgical equipment
- Stainless machines and equipment
- Rubber
- Plastic surfaces

Advantages

- Neutralizes alkaline and hard water residues
- Safe on stainless steel
- Prevents mineral buildup
- Prolongs instrument life
- Prevents spotting from alkaline cleaning
- Improves rinse water quality

Technical Data

- 2.1 – 3.1
- Bio-degradable
- Non-flammable
- Dilution: 1/4 to 1 oz product per 1 Gallon of water

Common Packaging

- 4 – 1 Gallon Bottles
- 2 – 2.5 Gallon Pails
- 5 Gallon Pails
- 33 Gallon Drums
- 55 Gallon Drums

